nutrifit / nutrisafe2 PVC, PUR & SIL feeding tubes

Ref nutrifit: 361 / 362 / 363 / 364 / 366 / 1361 / 1362 / 1363 / 1364 / 1366 / 2331 / 2332 / 2336 / 2395

Ref nutrisafe2: 361 / 362 / 363 / 364 / 1361 / 1362 / 1363 / 1364 / 1365 / 2331 / 2332 / 2395 / 80363 / 80362 / 81361 / 81363 / 81362 / 82332

English

PVC, PUR and SIL feeding tubes for safe enteral nutrition Nutrifit and Nutrisafe 2

Nutrifit and Nutrisafe 2 feeding tubes are indicated for enteral feeding only.

Centimeter numerical markings from 5 to 35 cm for the 40, 50 and 75 cm feeding tubes and markings from 5 to 70 cm for the 90, 125 and 160 cm feeding tubes.

Feeding tubes of the Nutrifit range are compliant with the ENFit secure connection from ISO 80369-3 requirements and must be used with dedicated disposable accessories.

Nutrisafe 2 feeding tubes are indicated for the enteral feeding in neonatology and pediatrics. They are fitted with the Nutrisafe 2 secure connection. They are compliant with ISO 80369-1 requirements and must be used with dedicated Nutrisafe 2 accessories.

Nutrifit and Nutrisafe 2 feeding tubes are available with a color coded obturator:

- Pre-mounted retaining strap color coded by size for Nutrifit.
- Pre-mounted obturator cap color coded by size for Nutrisafe2.

Available in PVC and polyurethane (PUR):

- Length 40, 50, 75, 90, 125 & 160 cm for Nutrifit
- Length 40, 50, 75, 125 et 160 cm for Nutrisafe 2

Available in silicone (SIL):

- Length 50, 90, 125 cm for Nutrifit
- Length 50 & 125 cm for Nutrisafe 2

It is recommended to replace Nutrifit / Nutrisafe 2 feeding tubes in PVC within 5 days and in PUR or Silicone within 4 weeks.

Tube placement

Gather all equipment required:

- Device to clean the nostrils.
 - Disposable gloves
 - Fixation device
 - If required a lubricant or stylet
 - Feeding tube or feeding tube with stylet (for size and material selection, refer to doctor prescriptions).
- Inform the patient (if this is possible).

Placement of the feeding tube

For oral or nasal placement, measure the length of tube to be inserted according to standard protocols.

Lubricate the feeding tube.

Introduce the feeding tube in one of the nostrils and move the tube carefully towards the posterior nasopharynx, oriented towards the oesophagus .

Proceed slowly to avoid the risk of the patient vomiting.

Bend the head forward to close the trachea and open the esophagus. Advance the tube into the esophagus.

Once the mark is reached, stop the advancement of the tube.

Fixation of the feeding tube:

After preparing the skin, dress adhesive strip around the tube and fix it on the nose.

Verification:

Check feeding tube positioning preferably by x-ray. If using a feeding tube with stylet, remove gently the stylet from the feeding tube.

Connection

The connection of Nutrifit feeding tubes is secured thanks to its ENFit connector from the ISO80369-3 norm.

The Nutrisafe 2 connection features specific design for Neonatal and paediatric use. It is compliant with ISO 80369-1.

So Nutrisafe2 helps to prevent misconnections with IV line (luer) and other ISO 80369 connectors (respiratory, urology, etc)

The connections are quick and simple. No change in technique is required, and specific training is not necessary. It is easy to connect the accessories, just screw the tab components.

Precaution for use

The rules of good hygiene should be followed at each step of the placement of the tube.

The placement of the tube is verified once a day and systematically before each use (in any case follow the protocols of your institution).

Before each use, check that the feeding tube adaptor is connected to an enteral port and not to an IV set.

Strictly follow the protocols of your institution to clean the tube between two uses
Nutrisafe 2 and Nutrifit feeding tubes are contraindicated for patients with congenital abnormalities of the gastrointestinal tract.

Feeding tubes less than 6 Fr are mainly used in neonatology and paediatrics

Feeding tubes from 8 Fr to 12 Fr are mainly used in paediatrics and for adult patients.

Feeding tubes of 14 Fr or more are mainly used for adult patients

In case of resistance, do not force the administration of the enteral nutrition formula.

For the feeding tubes with stylet:

Prior to insertion, always check the stylet does not come out through one of the eye of the tube. If it is the case, replace it correctly.
Never insert back the stylet in a feeding tube already placed in the patient.

Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

This device is not made with dry or natural rubber latex.

To dispose of, place the contaminated items in the appropriate receptacle.

This device incorporates connectors which are designed to reduce the likelihood of tubing misconnection. However, the potential to misconnect this device with medical devices of other healthcare applications still exists, including intravascular devices, neuraxial devices, hypodermic applications, breathing systems and driving gas devices, urinary devices and limb cuff inflation devices.

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